

A comparison of metagenome assembly strategies

Ino de Bruijn

Content

- Introduction on Assembly
 - Popular methods
- Benchmark Method
 - Which assemblers have been validated?
 - How have they been validated?
- Benchmark Results
 - What is the best metagenomic assembler?

Introduction on Assembly

- Different technologies to sequence your sample
 - Sanger (400-600nt, error rate 0.001-1.0%)
 - 454 (400nt, error rate 1-4%)
 - Illumina (100nt, error rate 0.1-1%)
 - PacBio (3k-5knt, error rate 13-20%)
 - Ion Torrent (200-400nt, error rate 0.5-2.5%)
- Assembly turns reads into contigs and/or scaffolds
- Algorithms fit sequencing technologies
- Three general approaches
 - Overlap-Layout-Consensus (Long reads Celera 2000, Newbler 2005)
 - Greedy (Short reads, SSAKE 2007)
 - de Bruijn Graph (Short reads, Velvet 2008)



Introduction on Assembly

- Choosing the right K...
 - larger K more specific less coverage (span repeats, regions occurring twice, less connections in the graph)
 - Smaller K more sensitive more coverage (more connections in the graph)
 - Ideally combine both
- De Bruijn Graph potential information available
 - Overlap between kmers
 - Kmer coverage (how often does a kmer occur)
 - Read that created the kmer (choose between paths)
 - Insert size distribution between pairs (if paired reads were used)
- Programs differ in
 - 1) how the graph is stored
 - 2) how the graph is traversed

Introduction on Assembly


- Metagenomics uses DNA from environmental sample
 - Not all microbes can be cultured in the lab
 - Study microbial communities in their natural environment
- Metagenomic assembly more difficult than single genome assembly
 - Number of genomes unknown (maybe a rough idea)
 - Coverage of genomes differs (different abundances of genomes)
 - Closely related strains complicate the graph (in de Bruijn: anything that shares stretches of DNA longer than K)

Benchmark Method

- Use an *in vitro* simulated metagenome with known species
 - 59 species, total size 195Mb
 - Two abundance distributions of genomes
 - Even, approximately the same genome copy numbers 
 - Uneven, log-normal distribution of phyla similar to soil 
 - Assemblies can be validated by aligning them against the reference genomes (we used nucmer for this)
- Use Illumina paired end reads since that is currently one of the most popular sequencing techniques for metagenomic samples

Benchmark Method

- Enormous amount of assembly strategies possible
 - Select a number of assembly strategies to test



Benchmark Method

- Enormous amount of assembly strategies possible
 - Select a number of assembly strategies to test

- Contiguing with




Velvet



Meta-Velvet



Ray



Benchmark Method

- Enormous amount of assembly strategies possible
 - Select a number of assembly strategies to test

- Contiguing with



Velvet



Meta-Velvet



Ray


- Merging with



Newbler



Minimus2



Benchmark Method

- Enormous amount of assembly strategies possible
 - Select a number of assembly strategies to test

- Contiguing with



Velvet



Meta-Velvet



Ray

- Merging with



Newbler



Minimus2

- Scaffolding with




Velvet



Meta-Velvet



Ray



Benchmark Method

- Enormous amount of assembly strategies possible
 - Select a number of assembly strategies to test

- Contiguing with



Velvet



Meta-Velvet



Ray

- Merging with



Newbler



Minimus2

- Scaffolding with



Velvet




Meta-Velvet



Ray




Bambus2




Benchmark Method

- Enormous amount of assembly strategies possible
 - Select a number of assembly strategies to test


- Contiging with



- Merging with



- Scaffolding with



Total of 21 strategies

Benchmark Method

- How to compare all these different strategies?
- Validation of metagenomic assembly often focuses on one or more of the following points:
 - 1) contig/scaffold length distribution
 - 2) contig/scaffold coverage of the reference metagenome
 - 3) chimericity and erroneoususness of the contigs/scaffolds
 - 4) functional annotation
 - 5) phylogenetic classification
- We focus on the first three since those also tend to improve 4 and 5 as shown by Mende et al (2012)
- Select winner in three categories


Benchmark Results

1) contig/scaffold length distribution

- Popular statistic: N50 length (or L50)
 - Weighted median of contig lengths (contigs weighted by length)
 - 50% of all bases in the assembly are in contigs \geq L50
 - Bigger is better
- Problem when comparing between different assemblers: cut-off length matters and sum of bases between assemblies differs
 - Advantage for assemblers that only output long contigs

Benchmark Results

1) contig/scaffold length distribution




Shows best single kmer strategy per assembler based on L50

Dots are L50 values at cut-off 100

Benchmark Results

1) contig/scaffold length distribution




Difficult to say which is winning based on length distribution alone.

Benchmark Results



1) contig/scaffold length distribution




Merging increases lengths and makes length distributions more similar.

L50 before → after

V 665 → 3487 / 2553


M-V 6892 → 43685 / 45085

R 919 → 4076 / 1731



Benchmark Results



1) contig/scaffold length distribution



Scaffolding most notable
with Velvet and Meta-
Velvet

L50 before → after

V	665	→	4453
M-V	6892	→	12604
R	919	→	1002





Benchmark Results

2) contig/scaffold coverage of the reference metagenome

- So the length distributions don't say a lot, how much of the reference metagenome is actually covered?

Benchmark Results

2) contig/scaffold coverage of the reference metagenome




Metagenome coverage:
Number of non-overlapping bases covered by the assembly expressed as ratio of entire metagenome

Each line is an assembly
Assemblies with highest L50 chosen

Benchmark Results


2) contig/scaffold coverage of the reference metagenome



Velvet does a good job covering the metagenome.


Benchmark Results

2) contig/scaffold coverage of the reference metagenome



Benchmark Results

2) contig/scaffold coverage of the reference metagenome



Merging increases length while keeping coverage high. Still original Velvet contigs win overall



Benchmark Results

3) chimericity and erroneoususness of the contigs/scaffolds

- If lengths increase, but coverage doesn't. Do we have duplicate contigs or erroneous contigs?

Benchmark Results

3) chimericity and erroneoususness of the contigs/scaffolds




$$\sum_{\text{contigs}} \frac{\# \text{ bases in best alignment contig}}{\text{total nr of bases assembly}}$$
Global purity:

Each dot is an assembly

Point size indicates kmer size (21-79)

Benchmark Results


3) chimericity and erroneousness of the contigs/scaffolds



Ergo, Merging and Bambus2 both perform well purity-wise






Benchmark Results

3) chimericity and erroneousness of the contigs/scaffolds



Ergo, Merging and Bambus2 both perform well purity-wise

Conclusions

- Difficult to name one true winner. Depends on what the post-processing will be like
 -  If short contigs useful: **Ray** or Velvet. Velvet preferred if really short contigs 100-1000 (more coverage of metagenome).
 -  If only interested in long contigs and less so in their purity, perhaps choose **Meta-Velvet**.
 -  Merging seems to work quite well although you lose some metagenome coverage
 - Might want to avoid scaffolding since it doesn't work that well
 -  bambus2 in particular loses a lot of metagenome coverage
 -  Internal scaffolding works better than bambus2, but merging alone does a better job

Future work

- Try different mock communities, preferably higher covered samples
- See how results hold with functional annotation
 - How many false positives and true positives versus reference metagenome
 - What is an adequate cut-off length to use for the assemblies?
- Compare reference-less validation against this validation
 - Compare for instance with Feature Response Curve (Vezzi, 2012)
- Make the validation pipeline available as a web service
 - Give reference and assembly, watch the performance of your assembly
 - Make your own assembly of existing reference and reads and compare against existing ones
 - Allow users to add their own analysis. For example through ipython notebook

Questions?